

# Welcome to AEcroscopy for Automated and Autonomous Microscopy

## Contents

- Introduction to AEcroscopy
- High Throughput Experimentation

Under construction (Here, we will introduce automated and autonomous experiment, the requirement of autonomous experiment, the objective of AEcroscopy... )

AEcroscopy is a cutting-edge Python package that integrates with our self-developed Application Programming Interface (API) BEPyAE, aiming to revolutionizing microscopy measurements. AEcroscopy enables easy implementation of experiments formulated with Python program.....allowing users to orchestrate complex experiments, design workflows, and gather data more efficiently.

In addition, AEcroscopy allows users to integrate machine learning models into the operating experiment, harnessing the power of AI for real-time data analysis and decision-making. This will enable rapid iterations in experiment and largely accelerate scientific discovery.

With AEcroscopy, a microscopy capable of automating operation, actively learning from on-the-fly results, adapting experiments, optimizing parameters, driving real-time decision-making... becomes a reality.

Discover the potential of AEcroscopy and experience the future of microscopy experimentation.

*Yongtao Liu, June 2023*

## Introduction to AEcroscopy

[Skip to main content](#)

Welcome to our comprehensive guide to the AEcroscopy for scanning probe microscopy (SPM). AEcroscopy is a dedicated Python package incorporating a wide range of commands and functions that facilitate basic microscope operations with BEPyAE.exe, enabling users to perform tasks with ease.

Here, we will introduce you to the essential Python commands and functions available within the AEcroscopy. We will explore their capabilities, syntax, and usage, empowering you to harness the full potential of the AEcroscopy-BEPyAE-SPM system. Whether you are a seasoned Python developer or new to the world of programming, experienced SPM users or new to the microscope society, this chapter will serve as a foundation for your journey into the exciting realm of Python-controlled SPM.

By utilizing AEcroscopy, you will gain the ability to write and execute SPM experiments in Python, revolutionizing the way we perform the microscope experiments.

*Yongtao Liu, June 2023*

## Guide to Essential Commands and Functionalities in AEcroscopy

*Yongtao Liu*

*June 2023*

Welcome to our comprehensive guide to the AEcroscopy for scanning probe microscopy (SPM). AEcroscopy is a dedicated Python package incorporating a wide range of commands and functions that facilitate basic microscope operations with BEPyAE.exe, enabling users to perform tasks with ease.

Here, we will introduce you to the essential Python commands and functions available within the AEcroscopy. We will explore their capabilities, syntax, and usage, empowering you to harness the full potential of the AEcroscopy-BEPyAE-SPM system. Whether you are a seasoned Python developer or new to the world of programming, experienced SPM users or new to the microscope society, this chapter will serve as a foundation for your journey into the exciting realm of Python-controlled SPM.

By utilizing AEcroscopy, you will gain the ability to write and execute SPM experiments in Python, revolutionizing the way we perform the microscope experiments.

### Install and Import

[Skip to main content](#)

First thing first, we need to install and import necessary packages, including AEcroscopy.

```
import os
import win32com.client
import numpy as np
import time
import h5py
import sidpy
import pyNSID
import matplotlib.pyplot as plt
from tqdm import tqdm

# import acquisition.py
from Acquisition_v0_9 import Acquisition # include the Acquisition_v.py in the same directory
```

## Start BEPyAE.exe and set VI

Then, we need to start the Python API BEPyAE.exe.

- Start BEPyAE.exe
- Set VI of BEPyAE; if this version includes PyScanner, also set VIs for PyScanner

```
newexp = Acquisition(exe_path = r"C:\Users\yla\Dropbox (ORNL)\My Files\AEcroscopy_BEPyAE\BEPyAE
```

## Initialize Igor AR18

Here, we connect BEPyAE with the microscope

- Set offline development. If you are doing offline development, set offline\_development=True, otherwise, if you are running microscope measurement, set offline\_development=False

Executing `init_BEPyAE()` command will also:

- Build a connection between BEPyAE and AR18
- Get parameters in AR18

```
newexp.init_BEPyAE(offline_development = True) # set offline_development=True if doing offline
                                                # executing this will also initialize AR18
```

[Skip to main content](#)


## Hereinafter

- if no parameter is input in a function, when you execute the function, it will take the default parameters in BEPyAE.exe
- Some functions print feedback after execution. This feedback can be turned off by setting feedbackon = False. You can turn off feedback when you include this function in some iterations.
- Note: Tip locations: -1 is the left handside for x-axis and bottom side for y-axis, 1 is the right handside for x-axis and top side for y-axis

## Set tip parameters

Here, we use *tip\_control()* function to set setpoint.

*tip\_control()* function also allows us to move tip to a specific location [next\_x\_pos\_00, next\_y\_pos\_01], e.g., in the below example, the tip will be moved to locaiton [0.5, 0.5]



```
newexp.tip_control(tip_parms_dict = {"set_point_V_00": 1, "next_x_pos_00": 0.5, "next_y_pos_01": 0.5,
                                      do_move_tip = True,
                                      do_set_setpoint = True,
                                      feedbackon=False) # Executing this code will set setpoint to 1 V,
                                         # and move tip to location [0.5, 0.5]
```

## Set IO

In order for the Python code to provide the correct commands, it is essential to provide the hardware components (e.g., AFM platform, voltage amplifiers, channel data, etc.) involved in our experiment. Thus, here we set information about hardwares with function *define\_io\_cluster()*.

```
newexp.define_io_cluster(IO_cluster_parms_dict = {"analog_output_amplifier_06": 1,  
                                         "channel_01_type_07": 1,  
                                         "channel_02_type_08": 2,  
                                         "channel_03_type_09": 3})
```

```
('0 Cypher AR18',  
 '6124',  
 4000000.0,  
 10.0,  
 10.0,  
 'AC and DC on A00',  
 1.0,  
 'none',  
 'none',  
 'none',  
 'external')
```

### BEPFM Measurement

## Set BE pulse parameters

One of the most popular experiment we can perform with AEcroscopy-BEPyAE.exe is the band excitation (BE) piezoresponse force microscopy experiments. Here we use *define\_be\_parms()* to set BE pulse parameters.

```
# set BE parameters  
newexp.define_be_parms(be_parms_dict = {"center_frequency_Hz_00": 335, "band_width_Hz_01": 100,  
                                         "amplitude_V_02": 1, "phase_variation_03": 1,  
                                         "repeats_04": 4, "req_pulse_duration_s_05": 4,  
                                         "auto_smooth_ring_06": 1},  
                         do_create_be_waveform = True, feedbackon=False)
```

After setting BE pulse parameters, we can start to do BE line measurement with the function `do_line_scan()`. In this function, we will need to provide the BE line scan pixel (i.e., `num_BE_Pulses_01`) and BE line scan locations (i.e. start location and stop location).

Note that

- This function is just a single BE line scan, not a raster image
- This function returns 5 datasets: quick\_fitting, complex spectra, and 3 channel data

```
# Do a single line scan
qk_fit, com_spec, chn1, chn2, chn3 = newexp.do_line_scan(line_scan_parms_dict = {"num_BE_pulses": 32,
                                                                 "start_x_pos_0": -0.5,
                                                                 "stop_x_pos_0": 0.5,
                                                                 "upload_to_daq": True, do_line_scan = T})
```

```
voltage offset and number of BE pulse are: (0.0, 32)
line scan start and end positions: (-0.5, 0.0, 0.5, 0.0)
```

## BE Raster Scan


We can also perform a square raster scan BE measurement using the function `raster_scan()`. In this function, we need to provide the raster scan pixel (i.e., `scan_pixel`) and scan region (i.e., scan start and stop points).

Note that

- `raster_scan` returns 3 sidpy datasets: BEPFM quick fitting, channels, and BE complex spectra
- `raster_scan` also saves these 3 sidpy dataset in a h5 file named `file_name`

```
# Do a 64*64 raster scan
dset_pfm, dset_chns, dset_cs = newexp.raster_scan(raster_parms_dict = {"scan_pixel": 32, "scan_x_start": -0.8, "scan_x_stop": 0.8, "scan_y_start": -0.8, "scan_y_stop": 0.8}, file_name = "raster.h5")

# if you see below error, check if you set IO channels manually---this has to be done manually
### TypeError: When specifying values over which a parameter is varied, values should not be a
```



```
[progress: 0:00:38] |***** | (ETA: 0:00:01) C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning: validate_h5_dimension may be removed in a future version',  
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning: validate_h5_dimension may be removed in a future version',  
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning: validate_h5_dimension may be removed in a future version',
```

## BE Raster Scan Results


We can visualize BE raster scan results including quick fitting images and channel images.

```
# Quick fit BE images
print(dset_pfm) # sidpy dataset of BE quick fit

# plot BEPFM quick fit data
f, (ax1, ax2, ax3, ax4) = plt.subplots(1, 4, figsize = (20, 5), dpi = 100)
ax1.imshow(dset_pfm[:, :, 0])
ax2.imshow(dset_pfm[:, :, 1])
ax3.imshow(dset_pfm[:, :, 2])
ax4.imshow(dset_pfm[:, :, 3])
```

```
sidpy.Dataset of type IMAGE_STACK with:  
dask.array<array, shape=(32, 32, 5), dtype=float64, chunksize=(32, 32, 5), chunktype=numpy.ndarray  
data contains: quick fit pfm (generic)  
and Dimensions:  
y axis: y axis (m) of size (32,)  
x axis: x axis (m) of size (32,)  
BE responses: channels (generic) of size (5,)
```


```
<matplotlib.image.AxesImage at 0x1ba0b485a50>
```



```
# Channel images  
print(dset_chns) # sidpy dataset of channels  
  
# plot channel data  
f, (ax1, ax2, ax3) = plt.subplots(1, 3, figsize = (15, 5), dpi = 100)  
ax1.imshow(dset_chns[0,:,:])  
ax2.imshow(dset_chns[1,:,:])  
ax3.imshow(dset_chns[2,:,:])
```


```
sidpy.Dataset of type IMAGE_STACK with:  
dask.array<array, shape=(3, 32, 32, 1), dtype=float64, chunksize=(3, 32, 32, 1), chunktype=nmpr  
data contains: channels (generic)  
and Dimensions:  
y axis: y axis (m) of size (3,)  
x axis: x axis (m) of size (32,)  
channels images: channels (generic) of size (32,)  
d: generic (generic) of size (1,)
```

```
<matplotlib.image.AxesImage at 0x1ba0b4448e0>
```



```
# sidpy dataset of complex spectra  
dset_cs
```

|              | Array        | Chunk         |
|--------------|--------------|---------------|
| <b>Bytes</b> | 0.97 MiB     | 0.97 MiB      |
| <b>Shape</b> | (32, 32, 62) | (32, 32, 62)  |
| <b>Count</b> | 1 Tasks      | 1 Chunks      |
| <b>Type</b>  | complex128   | numpy.ndarray |



We can also load the saved h5 file and analyze it after experiments

```
hf = h5py.File('test_0.hf5', 'r+')
sidpy.hdf_utils.print_tree(hf)
```


[Skip to main content](#)

```
/  
|   |- BE Channels  
-----  
|   |   |- Channels  
|   |   |   |- Channels  
|   |   |   |- channels images  
|   |   |   |- d  
|   |   |   |- x axis  
|   |   |   |- y axis  
|   |- BE Complex Spectra  
-----  
|   |   |- Complex Spectra  
|   |   |   |- Complex Spectra  
|   |   |   |- c  
|   |   |   |- location index x  
|   |   |   |- location index y  
|   |- BE Parameters  
-----  
|   |   |- frequency  
|   |   |- pulse parameters  
|   |   |- scan size  
|   |- BE Quick Fitting  
-----  
|   |   |- Quick Fitting  
-----  
|   |   |   |- BE responses  
|   |   |   |- Quick Fitting  
|   |   |   |- x axis  
|   |   |   |- y axis
```

```
be = hf["BE Quick Fitting/Quick Fitting/Quick Fitting"]
```

```
be = np.asarray(be)  
plt.imshow(be[:, :, 0])
```

```
<matplotlib.image.AxesImage at 0x1ba0c06efb0>
```



### BEPS Measurements

## Set BEPS parameters

In addition to BE piezoresponse force microscopy image measurement, we can also do BE piezoresponse spectroscopy (BEPS) measurement with the function `do_BEPS_xxx()`.

However, before doing BEPS measurement, we first need to define BEPS measurement parameters with the function `define_BEPS_parameters()`. In this function, we need to provide BEPS waveform parameters, such as amplitude, steps, cycle numbers, etc.

```
# Define BEPS parameters
newexp.define_BEPS_parameters(beps_parms_dict = {"amplitude_V_00": 8, "offset_V_01":0,
                                                 "read_voltage_V_02": 0, "step_per_cycle_03": 3,
                                                 "num_cycles_04": 3, "cycle_fraction_05": 0,
                                                 "cycle_phase_shift_06": 0, "measure_loops_07": 0,
                                                 "transition_time_s_08": 1E-3, "delay_after_step_09": 0})
```

[Skip to main content](#)

```
BEPS parameters are: (8.0, 0.0, 0.0, 32, 3, 0, 0, 0, 0.001, 0.0, 0.0, 0.0, 0.0, 1, 0, 8.0, 8.0, 7.0,
```

## Do grid BEPS

After defining BEPS parameters, we can execute `do_BEPS_grid()` function to perform a BEPS measurement. In the `do_BEPS_grid()` function, we can provide grid pixel number in `beps_grid_parms_dict()`


```
# do BEPS
beps_waveform, beps_quick_fit, beps_cx, beps_chns = newexp.do_beps_grid(beps_grid_parms_dict =
                                                               file_name = "BEPs_grid")
```

```
[progress: 0:00:39] |*****| (ETA: 0:00:01) C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
```

Plot the BEPS waveform used in the experiments.

```
f, ax = plt.subplots(figsize = (3, 2))
ax.plot(beps_waveform)
ax.set_ylabel("Voltage (V)")
ax.set_xlabel("Step")
```

```
Text(0.5, 0, 'Step')
```



[Skip to main content](#)

During experiments, BEPS loops can be visualized in *BEPyAE.exe*. After experiments, we can also plot all BEPS loops as below.

```
# convert sidpy dataset to numpy array
quick_fit = np.asarray(beps_quick_fit)

# if measure_loops_07 is on and off, we need to separate on field and off field responses.
on_field_quick_fit = np.zeros((quick_fit.shape[0], quick_fit.shape[1], len(beps_waveform)))
off_field_quick_fit = np.zeros((quick_fit.shape[0], quick_fit.shape[1], len(beps_waveform)))
for i in range (len(beps_waveform)):
    on_field_quick_fit[:, :, i] = quick_fit[:, :, 2*i]
    off_field_quick_fit[:, :, i] = quick_fit[:, :, 1+2*i]


##### Plot on field and off responses #####
idx = 0 # index of the loop to plot
f, (ax1, ax2, ax3, ax4) = plt.subplots(1, 4, figsize = (18, 3), dpi = 100)
f.subplots_adjust(wspace = 0.4)
ax1.set_title("Off Field, Amplitude vs. Voltage")
ax1.plot(beps_waveform, off_field_quick_fit[idx, 0, :])
ax1.set_xlabel("Voltage (V)")
ax1.set_ylabel("Amplitude (a.u.)")

ax2.set_title("Off Field, Phase vs. Voltage")
ax2.plot(beps_waveform, off_field_quick_fit[idx, 3, :])
ax2.set_xlabel("Voltage (V)")
ax2.set_ylabel("Phase (rad)")

ax3.set_title("On Field, Amplitude vs. Voltage")
ax3.plot(beps_waveform, on_field_quick_fit[idx, 0, :])
ax3.set_xlabel("Voltage (V)")
ax3.set_ylabel("Amplitude (a.u.)")

ax4.set_title("On Field, Phase vs. Voltage")
ax4.plot(beps_waveform, on_field_quick_fit[idx, 3, :])
ax4.set_xlabel("Voltage (V)")
ax4.set_ylabel("Phase (rad)")
```

Text(0, 0.5, 'Phase (rad)')



[Skip to main content](#)

## Do BEPS at a specific locations

In addition to BEPS measurement at grid locations, the `do_BEPS_specific()` function also enables BEPS measurement at specified locations by providing specified coordinates to `coordinates`

```
# do BEPS at specific locations
beps_waveform, beps_quick_fit, beps_cx, beps_chns = newexp.do_beps_specific(coordinates = np.asarray([100, 130]), file_name = "BEPS_a.h5")
```

However, when you want to measure a feature observed in an image, you often need to convert the location from the image to the coordinate for tip. Here we can use the `convert_coordinates()` function.

```
# e.g. convert the location [100, 130] at an image to the coordinate of tip.
locations = np.asarray([100, 130])
converted_coor = newexp.convert_coordinates(original_coordinates = locations, num_pix_x = 256, num_pix_y = 256)
print(converted_coor)
```

```
[-0.21875  0.015625]
```

After we have the coordinate for tip location, we can use the `do_BEPS_specific()` to perform BEPS measurement here

```
# do BEPS at specific locations
beps_waveform, beps_quick_fit, beps_cx, beps_chns = newexp.do_beps_specific(coordinates = converted_coor, file_name = "BEPS_a.h5")
```

```
[progress: 0:00:01] | (ETA: --::--) C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning: warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning: warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning: warn('validate_h5_dimension may be removed in a future version',
```

## Apply a pulse


[Skip to main content](#)

In addition to characterization, SPM also allows to manipulate objects. For example, applying a DC pulse via SPM tip to a ferroelectric film allows to flip the ferroelectric polarization.

In AEcroscopy, we can apply a pulse by using `define_apply_pulse()` function. We can set the pulse amplitude and pulse duration in this function, as shown below.

```
# Apply DC pulse
pulse, time = newexp.define_apply_pulse(pulse_parms_dict = {"pulse_init_amplitude_V_00": 0, "pu
                                         "pulse_final_amplitude_V_02": 0, "p
                                         "rise_time_s_05": 1E-3, "pulse_fina
                                         "pulse_repeats_06": 3},
                                         do_create_pulse = True, do_upload_pulse = True, do_appl
                                         # Plot applied pulse
f, ax = plt.subplots(figsize = (3, 2))
f.suptitle("Applied Pulse")
ax.plot(time, pulse)
ax.set_xlabel("Time (s)")
ax.set_ylabel("Voltage (V)")
```

Text(0, 0.5, 'Voltage (V)')



## Progress bar

AEcroscopy is designed for automated and autonomous microscopy measurements, so we also include a progress function `progress_bar()` that assists users to track the experiment progress. For example, when we perform 5 iterations BEPS measurement, we can use the progress bar to track the experiment progress as below.

```
iteration = 5
progress_bar = newexp.progress_bar(max_value = iteration)
for i in range (iteration):
    newexp.define_BEPS_parameters(do_VS_waveform = True, feedbackon = False)
    progress_bar.update(i)
```

```
[progress: 0:00:02] |*****| (ETA: 0:00:00)
```

### Capability with FPGA

In AECroscopy-BEPyAE.exe system, users have various options of scan trajectories when equipped with FPGA, such as spiral scan or any customized trajectories.

## Spiral Scan

Here we will show how to do spiral scan with the *fpga\_spiral\_scan()* function. Using this function, we set the spiral parameters, e.g., inner radius, outer radius, spiral cycles, duration, etc; the function will return the spiral result and save the result as a H5 file.

```
spiral_results = newexp.fpga_spiral_scan(spiral_parms_dict = {"spiral_inner_radius_x_V_00": 0,
                                                               "spiral_inner_radius_y_V_02": 0,
                                                               "spiral_N_cycles_04": 5, "spiral_"
                                                               "spiral_dose_distribution_06": 1,
                                                               "spiral_return_opt_08": 0, "scan_
                                                               "scan_y_offset_V": 0, "scan_rotat
do_scan_update = True, do_scan = True,
file_name = "spiral_scan")
```

## BE Spiral Scan

If we apply BE excitation waveform during spiral scan, we can also perform BE spiral measurement. The function to perform BE spiral measurement is *fpga\_spiral\_scan\_BE()*. The BE waveform parameters can be set either in this function or beforehand.

```
# set BE parameters, set spiral parameters, and do spiral BE
fpga_results, be_results = newexp.fpga_spiral_scan_BE(be_parms_dict = {"center_frequency_Hz_00": 335000.0, "amplitude_V_02": 1, "phi": 0.004, "repeats_04": 4, "req_pulses": 128, "do_scan_update": True, "do_BE_arb_line_scan_01": True, "spiral_recording": False}, spiral_parms_dict = {"spiral_inner_radius": 100000.0, "spiral_outer_radius": 100000.0, "spiral_inner_start": 0.0, "spiral_outer_start": 0.0, "spiral_N_cycles_04": 4}, num_BE_pulse = 128, do_scan_update = True, do_BE_arb_line_scan_01 = True, spiral_recording = False)
```

BE parameters are: (335000.0, 100000.0, 1.0, 1.0, 4, 0.004, 1, 3352.2952763920002, 0.1215945906)

## Tip Control by FPGA

FPGA can also drive the probe to a specific location with the function *fpga\_tip\_control()*. We need to input the tip location parameters when using this function, as shown below, executing the below command will move tip from [0, 0] to [0.5, 0.5].

```
newexp.fpga_tip_control(fpga_tip_parms_dict={"strat_x_position_V_00": 0, "strat_y_position_V_01": 0, "final_x_position_V_02": 0.5, "final_y_position_V_01": 0.5, "make_cur_pos_start_pos": False, "do_probe_move_update": True, "do_probe_move": True})
```

If we are moving the tip from current location, we can set the *make\_cur\_pos\_start\_pos = True* and provide the final location, as shown below.

```
newexp.fpga_tip_control(fpga_tip_parms_dict={"final_x_position_V_02": 0.5, "final_y_position_V_01": 0.5, "make_cur_pos_start_pos": True, "do_probe_move_update": True, "do_probe_move": True})
```

## FPGA Driven Line by Line Raster Scan

We can also use *fpga\_linebyline\_raster\_scan()* function to perform slow raster scan or line by line scan.

• when we set do\_full\_raster\_scan — Tip will wait to advance to next line — False and

[Skip to main content](#)

shown below.

```
raster_full = newexp.fpga_linebyline_raster_scan(line_by_line_raster_dict = {"raster_scan_size":  
    "raster_scan_size":  
    "raster_N_scan_lin":  
    "raster_line_durat:  
    "scan_x_offset_V_0":  
    "scan_y_offset_V_0":  
    "scan_rotation_deg":  
    initialize_line_by_line_raster=True, do_full_r:  
    wait_to_advance_to_next_line=False, do_next_ra:  
    stop_full_raster_scan=False})
```

- Alternatively, when we set *do\_full\_raster\_scan = False*, *wait\_to\_advance\_to\_next\_line = True*, and *do\_next\_raster\_line\_only = True*, this function can be used to perform scan line by line, as shown below.

```
raster_line = newexp.fpga_linebyline_raster_scan(line_by_line_raster_dict = {"raster_scan_size":  
    "raster_scan_size":  
    "raster_N_scan_lin":  
    "raster_line_durat:  
    "scan_x_offset_V_0":  
    "scan_y_offset_V_0":  
    "scan_rotation_deg":  
    initialize_line_by_line_raster=True, do_full_r:  
    wait_to_advance_to_next_line=True, do_next_ra:  
    stop_full_raster_scan=False})
```

# High Throughput Experimentation

This chapter is about high throughput experiment workflows

...Under construction (we will introduce high-throughput experiment here)...

*Yongtao Liu, June 2023*

## High Throughput Domain Writing Workflow

*Yongtao Liu*

- - - - -

[Skip to main content](#)

## Install and Import

```
import os
import win32com.client
import numpy as np
import time
import h5py
import sidpy
import pyNSID
import matplotlib.pyplot as plt
from tqdm import tqdm

# import acquisition.py
from Acquisition_v0_9 import Acquisition # include the Acquisition_v.py in the same directory
```

## Start BEPyAE.exe and set VI

- Start BEPyAE.ext
- Set VI of BEPyAE; if this version includes PyScanner, also set VIs for PyScanner

```
newexp = Acquisition(exe_path = r"C:\Users\yla\Dropbox (ORNL)\My Files\AEcroscopy_BEPyAE\BEPyAE"
# exe_path is the directory of BEPyAE;
```

## Initialize Igor AR18

- Set offline development
- Build a connection between BEPyAE and AR18
- Get parameters in AR18

```
newexp.init_BEPyAE(offline_development = True) # set offline_development=True if doing offline
# executing this will also initialize AR18
```

Hereinafter

- If no parameters are provided in a function, executing the function will utilize the default parameters within *BEPyAE.exe*.
- Certain functions provide feedback after execution which can be disabled by setting

[Skip to main content](#)


instances when it is not required.

- Note: For tip locations, -1 corresponds to the left-hand side on the x-axis and the bottom side on the y-axis, while 1 corresponds to the right-hand side on the x-axis and the top side on the y-axis.

## Set tip parameters

- set setpoint, tip locations

(-1, 1)                  (1, 1)



(-1, -1)                  (1, -1)

```
newexp.tip_control(tip_parms_dict = {"set_point_V_00": 1, "next_x_pos_00": -0.5, "next_y_pos_01": 0.5, "do_move_tip": True, "do_set_setpoint": True}) # Executing this code will set setpoint to 1 V, # and move tip to location [0.5, 0.5]
```

```
Setpoint is: 1.0  
Tip parameters are: (-0.5, 0.5, 0.2)
```

## Set IO

This defines IO parameters, such as AFM platform: AR18, amplifiers, channel data types, etc

```
newexp.define_io_cluster(IO_cluster_parms_dict = {"analog_output_amplifier_06": 1,  
                                              "channel_01_type_07": 1,  
                                              "channel_02_type_08": 2, "channel_03_type_09":
```

```
('0 Cypher AR18',  
 '6124',  
 4000000.0,  
 10.0,  
 10.0,  
 'AC and DC on A00',  
 10.0,  
 'topography',  
 'current',  
 'aux',  
 'external')
```

## Set BE pulse parameters

```
# set BE parameters  
newexp.define_be_parms(be_parms_dict = {"center_frequency_Hz_00": 335, "band_width_Hz_01": 100,  
                                         "amplitude_V_02": 1, "phase_variation_03": 1,  
                                         "repeats_04": 4, "req_pulse_duration_s_05": 4,  
                                         "auto_smooth_ring_06": 1},  
                         do_create_be_waveform = True)
```

BE parameters are: (335000.0, 100000.0, 1.0, 1.0, 4, 0.004, 1, 3352.2952763920002, 0.1215945906)

## Run a BE Line scan to test parameters

- This is a single BE line scan
- This returns 5 datasets: quick\_fitting, complex spectra, and 3 channels

```
# Do a single line scan  
qk_fit, com_spec, chn1, chn2, chn3 = newexp.do_line_scan(line_scan_parms_dict = {"num_BE_pulses":  
                                         "start_x_pos_0":  
                                         "stop_x_pos_02":  
                                         upload_to_daq = True, do_line_scan = T}
```

```
voltage offset and number of BE pulse are: (0.0, 32)
line scan start and end positions: (-0.5, 0.0, 0.5, 0.0)
```


Tests Done

## Experiment Starts

In this experiment, we begin by applying a DC pulse to switch the ferroelectric polarization. Subsequently, a BEPFM (Bias-Enhanced Piezoresponse Force Microscopy) measurement is conducted to image the domain structure.


1. To initiate the measurement process, we first need to determine the location for each individual measurement. There are two scenarios to consider:

- For each measurement, a new location is chosen, requiring a location array to record all the measurements as demonstrated below.




[Skip to main content](#)

- Alternatively, all measurements are conducted at the same location. In this case, the measurement location needs to be initialized each time a new measurement is started, as shown below.



- Prior to experiments, we also need to establish the DC pulse parameters including pulse magnitude  $A$  and pulse length  $t$ , as shown below. Again, there are two scenarios to consider here:



- The pulse parameters can be pre-defined, e.g., the parameter values can be uniformly distributed within a specified range or customized to suit the experimental requirements.
- The pulse parameters can be random values within a defined space, typically these random values is a uniform distribution across that space in principle.

## Experiment 1. Perform each measurement at a new location with pre-defined pulse parameters

Prior to experiment, set a directory to save data

```
os.chdir(r"C:\Users\yla\Dropbox (ORNL)\My Files\AEcroscopy_BEPyAE")
```

Step 1. Generate a location array

```
# All locations span across [start_point_x, end_point_x] in x-direction and [start_point_y, end_point_y] in y-direction
# There are num_x rows and num_y columns in the locations array

start_point_x = -0.9    # Define location array parameters
end_point_x = 0.9
start_point_y = -0.9
end_point_y = 0.9
num_x = 2
num_y = 2

# Generate location array
pos_x = np.linspace(-0.9, 0.9, num_x)
pos_y = np.linspace(-0.9, 0.9, num_y)
pulse_pos = np.meshgrid(pos_x, pos_y)
pulse_pos_x = pulse_pos[0].reshape(-1)
pulse_pos_y = pulse_pos[1].reshape(-1) # pulse_pos_x and pulse_pos_y are the coordinates of all locations

# Set BEPFM image size
img_size = 0.1

# Check
if img_size > np.abs(pos_x[0]-pos_x[1]):
    print ("Alert: there will be image overlap along x-direction")
elif img_size > np.abs(pos_y[0]-pos_y[1]):
    print ("Alert: there will be image overlap along y-direction")
else:
    print("{} locations are ready for experiments".format(len(pulse_pos_x)))
```

4 locations are ready for experiments

## Step 2. Establish pulse parameters

```
# uniformly distributed pulse parameters
min_voltage = 5
max_voltage = 9
Vdc_amp = np.linspace(min_voltage, max_voltage, num_x) # pulse magnitude

min_time_log = -4
max_time_log = 1
Vdc_time = np.linspace(min_time_log, max_time_log, num_y, dtype = np.float32())
Vdc_time = np.power(10, Vdc_time) # pulse time

# Establish pulse parameters
Vdc = np.meshgrid(Vdc_amp, Vdc_time)
Vdc_amp = Vdc[0].reshape(-1)
Vdc_time = Vdc[1].reshape(-1)

if len(Vdc_amp) > len(pulse_pos_x):
    print ("Error: No enough locations to test all pulse conditions")
else:
    print ("{} pulse parameters are ready for experiments".format(len(Vdc_amp)))

# save pulse condition
np.save("Vdc_list.npy", np.asarray([Vdc_amp, Vdc_time]))
```

4 pulse parameters are ready for experiments

### Step 3. Start experiment

```

for i in tqdm(range(len(Vdc_amp))):
    ##### Move tip to the pulse location #####
    newexp.tip_control(tip_parms_dict = {"set_point_V_00": 1,
                                         "next_x_pos_00": pulse_pos_x[i],
                                         "next_y_pos_01": pulse_pos_y[i]},
                        do_move_tip = True, do_set_setpoint = True)
    time.sleep(0.2)

    ##### Apply pulse #####
    V_amp = Vdc_amp[i]
    V_time = Vdc_time[i]
    newexp.define_apply_pulse(pulse_parms_dict = {"pulse_init_amplitude_V_00": 0, "pulse_mid_amplitude": 0, "pulse_final_amplitude_V_02": 0, "pulse_on_duration": 1E-4, "pulse_final_duration": 1, "pulse_repeats_06": 1},
                               do_create_pulse = True, do_upload_pulse = True, do_apply_pulse =
    #
    time.sleep(1)
    newexp.define_apply_pulse(pulse_parms_dict = {"pulse_init_amplitude_V_00": 0, "pulse_mid_amplitude": 0, "pulse_final_amplitude_V_02": 0, "pulse_on_duration": 1E-4, "pulse_final_duration": 1, "pulse_repeats_06": 1},
                               do_create_pulse = True, do_upload_pulse = True, do_apply_pulse =
    time.sleep(2)

    ##### Do BEPFM to image domain #####
    dset_pfm, dset_chns, dset_cs = newexp.raster_scan(raster_parms_dict = {"scan_pixel": 16,
                                                                           "scan_x_start": pulse,
                                                                           "scan_y_start": pulse,
                                                                           "scan_x_stop": pulse,
                                                                           "scan_y_stop": pulse,
                                                                           "file_name": "Domain_Writing_{}".format(i),
                                                                           "progress_on": False, "ploton": False})
    time.sleep(0.5)


# Plot BEPFM images
f, (ax1, ax2, ax3, ax4, ax5, ax6) = plt.subplots(1, 6, figsize = (30, 5), dpi = 100)
ax1.imshow(dset_pfm[:, :, 0])
ax2.imshow(dset_pfm[:, :, 1])
ax3.imshow(dset_pfm[:, :, 2])
ax4.imshow(dset_pfm[:, :, 3])
ax5.imshow(dset_chns[0, :, :])
ax6.imshow(dset_chns[1, :, :])
plt.show()

```

0% |

Setpoint is: 1.0  
Tip parameters are: (-0.9, -0.9, 0.5)

```
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
```




25% |

| 1/4 |

Setpoint is: 1.0  
Tip parameters are: (0.9, -0.9, 0.5)

```
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
```




50% |

| 2/4 |

[Skip to main content](#)

Tip parameters are: (-0.9, 0.9, 0.5)

```
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
```




75%

| 3/4 |

Setpoint is: 1.0

Tip parameters are: (0.9, 0.9, 0.5)

```
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
```



100%

| 4/4 |

Step 4. Do a BEPFM at the whole experiment area

[Skip to main content](#)


```

dset_pfm, dset_chns, dset_cs = newexp.raster_scan(raster_parms_dict = {"scan_pixel": 32, "scan_x_start": -1.0, "scan_x_stop": 1.0, "scan_y_start": -1.0, "scan_y_stop": 1.0}, file_name = "pfm_whole", ploton = False)

f, (ax1, ax2, ax3, ax4, ax5, ax6) = plt.subplots(1, 6, figsize = (30, 5), dpi = 100)
ax1.imshow(dset_pfm[:, :, 0])
ax2.imshow(dset_pfm[:, :, 1])
ax3.imshow(dset_pfm[:, :, 2])
ax4.imshow(dset_pfm[:, :, 3])
ax5.imshow(dset_chns[0, :, :])
ax6.imshow(dset_chns[1, :, :])
plt.show()

```

[progress: 0:00:35] |\*\*\*\*\* | (ETA: 0:00:01) C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf\_utils.py:376: FutureWarning: validate\_h5\_dimension may be removed in a future version', C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf\_utils.py:376: FutureWarning: validate\_h5\_dimension may be removed in a future version', C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf\_utils.py:376: FutureWarning: validate\_h5\_dimension may be removed in a future version',



## Experiment 2. Perform measurement with random pulse parameters

Prior to experiment, set a directory to save data

```
os.chdir(r"C:\Users\yla\Dropbox (ORNL)\My Files\AEcroscopy_BEPyAE")
```

Step 1. Generate a location array

[Skip to main content](#)

```

# All locations span across [start_point_x, end_point_x] in x-direction and [start_point_y, end_
# There are num_x rows and num_y columns in the locations array

start_point_x = -0.9    # Define location array parameters
end_point_x = 0.9
start_point_y = -0.9
end_point_y = 0.9
num_x = 2
num_y = 2

# Generate location array
pos_x = np.linspace(-0.9, 0.9, num_x)
pos_y = np.linspace(-0.9, 0.9, num_y)
pulse_pos = np.meshgrid(pos_x, pos_y)
pulse_pos_x = pulse_pos[0].reshape(-1)
pulse_pos_y = pulse_pos[1].reshape(-1) # pulse_pos_x and pulse_pos_y are the coordinates of all locations

# Set BEPFM image size
img_size = 0.1

# Check
if img_size > np.abs(pos_x[0]-pos_x[1]):
    print ("Alert: there will be image overlap along x-direction")
elif img_size > np.abs(pos_y[0]-pos_y[1]):
    print ("Alert: there will be image overlap along y-direction")
else:
    print("{} locations are ready for experiments".format(len(pulse_pos_x)))

```

4 locations are ready for experiments

## Step 2. Establish pulse space

```

# pulse magnitude space
min_voltage = 5
max_voltage = 9

# pulse time space
min_time_log = -4
max_time_log = 1

```

## Step 3. Start experiment

[Skip to main content](#)

```

# Set measurement iteration number
iterations = 100
# Set lists to save pulse parameters
Vdc_amp = []
Vdc_time = []

# Measurement starts
for i in tqdm(range(iterations)):
    ##### Move tip to the pulse location #####
    newexp.tip_control(tip_parms_dict = {"set_point_V_00": 1,
                                          "next_x_pos_00": pulse_pos_x[i],
                                          "next_y_pos_01": pulse_pos_y[i]},
                        do_move_tip = True, do_set_setpoint = True)
    time.sleep(0.2)

    ##### Apply pulse #####
    # Random pulse parameters
    V_amp = np.random.uniform(min_voltage, max_voltage)
    V_time = np.random.uniform(min_time_log, max_time_log)
    V_time = np.power(10, V_time)
    # Add parameters to Vdc lists
    Vdc_amp.append(V_amp)
    Vdc_time.append(V_time)

    # Apply pulse
    newexp.define_apply_pulse(pulse_parms_dict = {"pulse_init_amplitude_V_00": 0, "pulse_mid_amplitude_V_01": 0, "pulse_final_amplitude_V_02": 0, "pulse_on_duration_s_03": 1E-4, "pulse_final_duration_s_04": 1, "pulse_repeats_05": 1},
                               do_create_pulse = True, do_upload_pulse = True, do_apply_pulse =
                               #
                               time.sleep(1)
    newexp.define_apply_pulse(pulse_parms_dict = {"pulse_init_amplitude_V_00": 0, "pulse_mid_amplitude_V_01": 0, "pulse_final_amplitude_V_02": 0, "pulse_on_duration_s_03": 1E-4, "pulse_final_duration_s_04": 1, "pulse_repeats_05": 1},
                               do_create_pulse = True, do_upload_pulse = True, do_apply_pulse =
                               #
                               time.sleep(2)

    ##### Do BEPFM to image domain #####
    dset_pfm, dset_chns, dset_cs = newexp.raster_scan(raster_parms_dict = {"scan_pixel": 16,
                                                                           "scan_x_start": pulse_pos_x[i],
                                                                           "scan_y_start": pulse_pos_y[i],
                                                                           "scan_x_stop": pulse_pos_x[i+1],
                                                                           "scan_y_stop": pulse_pos_y[i+1],
                                                                           "file_name": "Domain_Writing_{}".format(i),
                                                                           "progress_on": False, "ploton": False})
    time.sleep(0.5)

    # Plot BEPFM images
    f, (ax1, ax2, ax3, ax4, ax5, ax6) = plt.subplots(1, 6, figsize = (30, 5), dpi = 100)
    ax1.imshow(dset_pfm[:, :, 0])

```

[Skip to main content](#)


```
ax4.imshow(dset_pfm[:, :, 3])
ax5.imshow(dset_chns[0, :, :])
ax6.imshow(dset_chns[1, :, :])
plt.show()

# Save pulse parameters
np.save("Vdc_list.npy", np.asarray([Vdc_amp, Vdc_time]))
```

0%|

Setpoint is: 1.0  
Tip parameters are: (-0.9, -0.9, 0.5)

```
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
```




1%|■

| 1/100 [e]

Setpoint is: 1.0  
Tip parameters are: (0.9, -0.9, 0.5)

```
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
```




2%|■

| 2/100 [01

[Skip to main content](#)

Tip parameters are: (-0.9, 0.9, 0.5)

```
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
```




3% |

| 3/100

Setpoint is: 1.0

Tip parameters are: (0.9, 0.9, 0.5)

```
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
C:\Users\yla\AppData\Local\anaconda3\lib\site-packages\pyNSID\io\hdf_utils.py:376: FutureWarning
  warn('validate_h5_dimension may be removed in a future version',
```



4% |

| 4/100

IndexError

```
Cell In[19], line 11
    7 # Measurement starts
    8 for i in tqdm(range(iterations)):
    9     . . . . .
```

Traceback (most recent call last)

[Skip to main content](#)

```
--> 11             "next_x_pos_00": pulse_pos_x[i],  
12             "next_y_pos_01": pulse_pos_y[i]},  
13         do_move_tip = True, do_set_setpoint = True)  
14     time.sleep(0.2)  
15 #####----- Apply pulse -----#####  
16 # Random pulse parameters
```

IndexError: index 4 is out of bounds for axis 0 with size 4

## Step 4. Do a BEPFM at the whole experiment area

```
dset_pfm, dset_chns, dset_cs = newexp.raster_scan(raster_parms_dict = {"scan_pixel": 32, "scan_  
"scan_y_start": -1.0, "sc  
"scan_y_stop": 1.0},  
file_name = "pfm_whole", ploton = False)  
  
f, (ax1, ax2, ax3, ax4, ax5, ax6) = plt.subplots(1, 6, figsize = (30, 5), dpi = 100)  
ax1.imshow(dset_pfm[:, :, 0])  
ax2.imshow(dset_pfm[:, :, 1])  
ax3.imshow(dset_pfm[:, :, 2])  
ax4.imshow(dset_pfm[:, :, 3])  
ax5.imshow(dset_chns[0, :, :])  
ax6.imshow(dset_chns[1, :, :])  
plt.show()
```

## Notebook for high throughput domain writing analysis

*YongtaoLiu*

*June15, 2023*

```
# You may need to sidpy if you did not do before  
# !pip install sidpy
```

## Import

[Skip to main content](#)

```
import numpy as np
import matplotlib.pyplot as plt
import os
import h5py
import sidpy
import cv2
import imutils
```

## Set the directory and root name of your dataset

```
# directory
path = r"C:\Users\yla\Dropbox (ORNL)\My Files\AEcroscopy_BEPyAE\AEcroscopy_data\Experiment1"

# name
file_name = "Domain_Writing_"

# spot numbers
num_x = 8
num_y = 8
```

## Load all data

```

count_imgs = 64 # variable of how many images you have in the directory
pixel = 64 # pixel of your image

# change working directory
os.chdir(path)

# create arrays for all data
amplitude = np.zeros((count_imgs, pixel, pixel))
phase = np.zeros((count_imgs, pixel, pixel))
frequency = np.zeros((count_imgs, pixel, pixel))
qfactor = np.zeros((count_imgs, pixel, pixel))
topography = np.zeros((count_imgs, pixel, pixel))

for i in range (count_imgs):
    h5 = h5py.File('Domain_Writing_{}_{0}.hf5'.format(i), 'r+')
    be_qf = h5["BE Quick Fitting/Quick Fitting/Quick Fitting"]
    be_ch = h5["BE Channels/Channels/Channels"]
    amplitude[i,] = be_qf[:, :, 0]
    phase[i,] = be_qf[:, :, 3]
    frequency[i,] = be_qf[:, :, 1]
    qfactor[i,] = be_qf[:, :, 2]
    topography[i,] = be_ch[0, :, :, 0]

nor_amplitude = (amplitude - amplitude.min()) / amplitude.ptp()
nor_phase = (phase - phase.min()) / phase.ptp()
nor_frequency = (frequency - frequency.min()) / frequency.ptp()
nor_topography = (topography - topography.min()) / topography.ptp()

```

## Plot data

Plot all results

```
for i in range (count_imgs):

    fig, axs = plt.subplots(1, 4, figsize=(16, 4), dpi = 100)
    fig.subplots_adjust(left=0.02, bottom=0.06, right=0.95, top=0.99, wspace=0.25)
    cm = 'viridis'
    shrink = 0.6


    im0 = axs[0].imshow(amplitude[i, ]*1000, origin = "lower", interpolation='nearest', cmap=cm)
    fig.colorbar(im0, ax=axs[0], shrink = shrink, label = "Amplitude (a.u.)")
    axs[0].axis("off")

    im1 = axs[1].imshow(phase[i, ], origin = "lower", interpolation='nearest', cmap=cm)
    fig.colorbar(im1, ax=axs[1], shrink = shrink, label = "Phase (rad)")
    axs[1].axis("off")


    im2 = axs[2].imshow(frequency[i, ]/1000, origin = "lower", interpolation='nearest', cmap=cm)
    fig.colorbar(im2, ax=axs[2], shrink = shrink, label = "Frequency (kHz)")
    axs[2].axis("off")

    im3 = axs[3].imshow(qfactor[i, ], origin = "lower", vmin = 0, vmax = 250, interpolation='nearest')
    fig.colorbar(im0, ax=axs[3], shrink = shrink, label = "Q Factor")
    axs[3].axis("off")


plt.show()
```




[Skip to main content](#)




[Skip to main content](#)




[Skip to main content](#)




[Skip to main content](#)




[Skip to main content](#)




[Skip to main content](#)




[Skip to main content](#)



[Skip to main content](#)



[Skip to main content](#)



Plot a specific channel together

```
fig, axes = plt.subplots(num_y, num_x, figsize=(num_x, num_y),
                       subplot_kw={'xticks':[], 'yticks':[]},
                       gridspec_kw=dict(hspace=0.02, wspace=0.02))

for ax, i in zip(axes.flat, range(count_imgs)):
    ax.imshow(amplitude[i], origin = "lower") # We are plotting amplitude now, you can change
```

